

Recording Concrete Use Case(UC3):

Actors:Person using MCT device

Precondition:MCT device is turned on

Description:

- 1)User selects programs or frequency by pressing the Ok button
- 2)Then chooses one of the 4 therapies or the 4 frequencies that are available by scrolling with the up and down arrow buttons
- 3)Once a user has selected a therapy or frequency a window is then opened asking the user if they would like to start the recording
- 4)The user selects the option to start recording the therapy or frequency
- 5)Once a user has finished recording the user ends the recording by pressing the back arrow button or menu button to leave the page.
- 6)Once finished recording the record is added to the history of treatment

Post condition:User has recorded a therapy and has added the recording to the history of treatment