

### Recording Concrete Use Case(UC3):

Actors:Person using MCT device

Precondition:MCT device is turned on

#### Description:

1)User goes to the tab bar and toggles the record to be turned on so that it starts recording

2)User then selects programs by pressing the Ok button

3)Then chooses one of the 4 therapies available by scrolling with the up and down arrow buttons

4)Once user has finished with therapy they go to the tab bar and toggle the record to be turned off so that the recording has finished

5)Once finished recording the record is added to the history of treatment

Post condition:User has recorded a therapy and has added the recording to the history of treatment